

Improving Outcomes In Cranioplasty – Clinical Results From 1055 Cranioplasties Using OssDsign® Cranial PSI

Abstract

Reconstruction of cranial defects can be a complex surgical procedure associated with an underestimated morbidity. This report describes the outcome of 1055 cranioplasties using OssDsign Cranial PSI, a patient-specific implant made from a calcium phosphate material reinforced with 3D printed titanium. All data was collected as part of post-market surveillance following introduction of the product in Europe, US and selected Asian markets in compliance with MEDDEV 2.7/1 rev.4 and MDR 2017/745. At a median follow up time of 21 months (range, 0-70 months) 44 implants (4.2%) had been explanted; whereof 22 (2.1%) of the implants were explanted due to early postoperative infections, 13 (1.2%) due to persistent wound dehiscence, 1 (0.1%) due to unsatisfactory aesthetical outcome and 8 (0.8%) due to other reasons, such as tumor recurrence.

Histological analysis of several implants explanted \geq 9 months following surgery revealed bony integration between the implant and the native bone, as well as new bone formation within and around the remaining calcium phosphate material.

Introduction

Cranioplasty is sometimes perceived as a straightforward procedure but literature confirms the opposite. Reconstruction of the cranium, especially involving large cranial defects, has high morbidity, regardless of the choice of reconstructive material. Autologous bone flaps have been the gold standard for a long time but rates of bone resorption and infection are high^{1,2}. Use of inert alloplastic materials such as titanium, PEEK or PMMA, tailored to the patients defect anatomy can be used but these materials may be less than optimal as they seem to be linked to high rates of implant exposure, infection and ultimately implant removal^{3,4}. Known risk factors of implant failure include irradiation, previous cranioplasty failures, thin and fragile soft tissue, exposed sinus cavities, age, and previous infections⁵.

OssDsign Cranial PSI is the only patient-specific cranial implant that combines mechanical performance with long-term bone integration and remodelling. The implant consists of a 3D printed medical grade titanium mesh skeleton, encased in a calcium phosphate material with clinical and pre-clinical evidence of bone regenerative characteristics. OssDsign Cranial PSI is designed to be used for non-load bearing applications in patients where cranial growth is complete, and for use with an intact dura with or without duraplasty. The device is custom-made to fit each patient-specific cranial defect.

To date (September 2020), a total of 1069 devices have been delivered to 110 European hospitals, 11 Asian hospitals and 60 hospitals in the USA.

Materials and methods

Post-market surveillance data collection as part of regulatory requirements has been continuously performed by OssDsign in compliance with MEDDEV 2.7/1 rev.4 and MDR 2017/745. Specific patient-related information, such as age, sex and underlying pathology is not applicable, as this is not revealed during the normal implant ordering process. This data presents the outcome of 1055 cases of cranial reconstruction using OssDsign Cranial PSI. 14 of the 1069 devices originally ordered

TABLE 1.

Device size (cm ²)	Number of devices (%)
\leq 50	180 (17)
51-100	225 (21)
101-150	418 (39)
151-200	219 (20)
201-250*	127 (3)
<i>Total</i>	<i>1055 (100)</i>

Size distribution of OssDsign Cranial PSI in clinical use.

*Devices shipped to Europe and Israel.

were not implanted for patient-specific reasons, none of which were device-related.

Of the 1069 OssDsign Cranial devices originally delivered, 76.5% (818/1069) were ordered by university hospitals with a high-level trauma unit (Table 2).

TABLE 2.

Hospital type	Number of hospitals
University Hospital	109
General Hospital	63
Army Hospital	1
Veteran Hospital	1
Private Hospital	7
<i>Total</i>	<i>181</i>

Hospital systems using OssDsign Cranial to date (September 2020).

One patient experienced a tumor recurrence 31 months following reconstructive surgery. This allowed for explantation of the implant and subsequent preparation of histological samples for analysis of bone formation.

Results

As per September 2020, a total of 1055 OssDsign Cranial PSI devices had been successfully implanted in US, European and Asian patients. At a median follow up time of 21 months (range 0-70 months) 44 implants (4.2%) had been explanted; whereof 22 (2.1%) of the implants were explanted due to early postoperative infections, 3 (0.3%) due to tumor recurrence, 13 (1.2%) due to persistent wound dehiscence, 3 (0.3%) due to early postoperative hematomas, 1 (0.1%) due to a non-device related dura rift post-op, 1 (0.1%) due to progression of autologous bone flap resorption and 1 (0.1%) due to unsatisfactory aesthetical outcome. (Table 3). None of the explantations were performed due to complications that were determined to be device related by the operating surgeon. Peer-reviewed data from a subpopulation of the material has shown similar figures in terms of device explantation due to infection⁶.

A handful of cases have been reported where the mosaic structure of the OssDsign Cranial PSI is visible through the skin. These patients have in general thin and fragile soft tissue. Even though no complications have occurred in these patients the unsatisfactory esthetical outcome should be reviewed, and potential design mitigations considered.

TABLE 3.

Primary cause of explantation	Number of patients (%)
Infection (Early post-op)	22 (2.1)
Tumor Recurrence	3 (0.3)
Persistent Wound Dehiscence	13 (1.2)
Hematoma (Early post-op)	3 (0.3)
Non-device Related Dura Rift	1 (0.1)
ABF Resorption	1 (0.1)
Aesthetic	1 (0.1)

Reasons for explantation of OssDsign Cranial PSI.

Histological analysis of one retrieved implant (Figure 1-2) showed that the calcium phosphate was partly transformed into new, well-vascularised osteonal bone after 31 months, indicating that the triphasic calcium phosphate composition has osteoconductive properties and that new bone growth can bridge between the ceramic tiles. This is consistent with earlier published data on use of the exact same calcium phosphate composition for cranial reconstruction^{6,7,8}. The regenerative features of the material has also been confirmed in several preclinical studies.




Figure 1. Histological evidence of bone formation at 31 months post implantation. Paragon stained sectioning of OssDsign Cranial PSI shows bony integration between the implant and the recipient bone (rb) as well as new bone formation (nb) within, and around the remaining calcium phosphate material (CaP) and supporting titanium structure (Ti).




Figure 2. Magnification of interface between recipient bone and OssDsign Cranial PSI following 31 months of implantation. The magnified picture clearly shows the viable new bone (nb) growing in and around the remnants of the triphasic calcium phosphate material (CaP) of OssDsign Cranial PSI. The interface between new bone and recipient bone (rb) shows complete integration of the implant.

A 52-week preclinical implantation study in a sheep model revealed the same pattern of host bone integration of the implant along with new bone formation in and around the calcium phosphate material^{9,10}.

In conclusion, OssDsign Cranial PSI has shown exceptional performance with an infection rate warranting implant removal of 2.1% in a patient population having undergone 1055 cranioplasties at a median follow up time of 21 months

The bone-regenerative capacity of the calcium phosphate material has been substantiated in preclinical studies and is supported by clinical experiences in multiple cases.

References

1. Korhonen TK, Salokorpi N, Niinimäki J, Serlo W, Lehenkari P, Tetri S. Quantitative and qualitative analysis of bone flap resorption in patients undergoing cranioplasty after decompressive craniectomy. *J Neurosurg*. 2018;1-10. doi:10.3171/2017.8.JNS171857.
2. Schwarz F, Dünisch P, Walter J, Sakr Y, Kalff R, Ewald C. Cranioplasty after decompressive craniectomy: is there a rationale for an initial artificial bone-substitute implant? A single-center experience after 631 procedures. *J Neurosurg*. 2016;124(March):710-715. doi:10.3171/2015.4.JNS159.710.
3. Thien A, King NKK, Ang BT, Wang E, Ng I. Comparison of polyetheretherketone and titanium cranioplasty after decompressive craniectomy. *World Neurosurg*. 2015;176-180. doi:10.1016/j.wneu.2014.06.003.
4. Coulter IC, Pesic-Smith JD, Cato-Addison WB, et al. Routine but risky: A multi-centre analysis of the outcomes of cranioplasty in the North-east of England. In: *Acta Neurochirurgica*. ; 2014.

5. Punchak M, Chung LK, Lagman C, et al. Outcomes following polyetheretherketone (PEEK) cranioplasty : Systematic review and meta-analysis. J Clin Neurosci. 2017;41:30-35. doi:10.1016/j.jocn.2017.03.028.
6. Kihlström Burenstam Linder K, Birgersson U, Lundgren K, et al. Patient-specific Titanium-Reinforced Calcium Phosphate Implant for the Repair and Healing of Complex Cranial Defects. World Neurosurg. 2019 Feb; 122:e399-e407. doi:10.1016/j.wneu.2018.10.061. Epub 2018 Oct 18.
7. Engstrand T, Kihlström L, Neovius E, et al. Development of a bioactive implant for repair and potential healing of cranial defects. J Neurosurg. 2014;120(1):273-277. doi:10.3171/2013.6.JNS1360.
8. Engstrand T, Kihlström L, Lundgren K, Trobos M, Engqvist H, Thomsen P. Bioceramic Implant Induces Bone Healing of Cranial Defects. Plast Reconstr surgery Glob open. 2015;3(8):e491. doi:10.1097/GOX.0000000000000467.
9. Gallinetti S, Burenstam Linder LK, Åberg J, et al. Titanium reinforced-calcium phosphate improves bone formation and osteointegration in ovine calvaria defects: a comparative 52-weeks study. Biomed Mater. 2020. doi: 10.1088/1748-605X/abca12.
10. Omar O, Engstrand T, Kihlström Burenstam Linder L, et al. In situ bone regeneration of large cranial defects using synthetic ceramic implants with a tailored composition and design. PNAS. 2020 117 (43): 26660-26671; doi: 10.1073/pnas.2007635117

About OssDsign Cranial PSI

OssDsign Cranial is a patient-specific implant based on a biocompatible calcium phosphate composition with a strong titanium skeleton embedded in the core of its ceramic tiles.

OssDsign Cranial PSI is intended for the reconstruction of cranial defects. It is indicated for non-load bearing applications for patients in whom cranial growth is complete, and for use with an intact dura, with or without duraplasty. Always read instructions for use which accompany the product for indications, contraindications, warnings and precautions.

About OssDsign

OssDsign is an innovator, designer and manufacturer of personalized bone replacement technology for cranial repair. We are surgeons, scientist and engineers – committed to improving outcomes in cranioplasty. For more information visit ossdsign.com.

