

OssDsign® Cranial PSI

PATIENT INFORMATION

GENERAL INFORMATION

OssDsign Cranial PSI is a medical device for reconstruction of cranial bone defects. Cranial PSI is intended for permanent implantation. The device consists of a 3D printed titanium mesh that is largely covered by a calcium phosphate ceramic. The ceramic part is approximately ~90% of the implant volume.

PRECAUTIONS

Magnetic Resonance Imaging (MRI)

OssDsign Cranial PSI is MR conditional. If you are to be subjected to magnetic resonance imaging, inform your healthcare

provider that you have an OssDsign Cranial PSI (patient-specific implant). Detailed information can be found here: www.ossdsign.com/OSSDSIGN_Cranial.

Radiotherapy

OssDsign Cranial PSI may affect radiotherapy dosage. If you are to be subjected to radiotherapy, inform your healthcare provider that you have an OssDsign Cranial PSI (patient-specific implant).



Always read the Instructions for Use which accompany the product for indications, contraindications, warnings and precautions.